

---

# **Doak MK4 Portable Surgical Table Features and Accessories**

---

# Doak MK4 Portable Surgical Table

---



## Doak Table accessories list

---



### Standard accessories:

- 2 fully adjustable arm boards
- 2 telescoping stainless steel IV poles
- 1 stainless steel Mayo tray
- 1 side-mounted surgical instrument tray
- 1 fluid containment sheet
- 1 folding lower shelf
- 2 anesthetist's equipment brackets
- 1 end-mounted anesthetist's equipment tray
- Soft carry/transport cases

### Optional accessories:

- Pelican 1780 and 1730 cases
- Run-flat, all-terrain wheels with casters
- Tempus Pro vital signs mounting bracket

## Doak Table with litter in Trendelenburg position

---



Trendelenburg or reverse Trendelenburg position offered from either end.


## Litter in Trendelenburg position – 15" height elevation

---



## Litter yoke and polypropylene retention straps

---





## Accepts Talon II and STANAG 2040 litters

---



**Folded dimensions: 41" x 15" x 12"**


---





## Stainless steel tie-downs for securing table for use during transportation

---



## Height-adjustable legs with self-leveling, non-skid feet

---



Operating height adjustable  
from 28" to 36"

**Fabricated with anodized, 6061 aluminum alloy tubing**

---



## North American or EU side rail with accessory clamp


---






## Fully adjustable arm boards and anesthetist's tray / equipment brackets

---



## Side-mounted surgical instrument tray


---





## Snap-on fluid containment sheet and folding lower storage shelf

---



## Swiveling, height-adjustable Mayo tray

---



**7.5 cm x 25 cm run-flat, multi-surface, locking wheels with custom casters**

---



**NOTE: Doak MK4 comes standard with four (4) non-skid, self-leveling feet. Wheels available as an option.**

## Doak Table with accessories prepared for transportation/storage

---





## Padded Soft Case with Shoulder Straps



*NOTE: Photo is of pre-production prototype*

## Optional Pelican 1780 and 1730 cases

---

