

## Features:

- Complete line of Flowmeters from Neonatal to Adult
- Available in Standard Thorpe style and Dial Style
- Standard 1/8" FNPT connection, but may also be fitted with different adapters and accessories to customize your needs
- MRI options are available
- Available in USA and ISO colors
- Smallest and lightest chromed brass bodied Flowmeter in the industry
- Color coded inner tube and labeling
- Easy to service
- No strip "D" shape knob
- Includes an inlet filter and inlet extension

**5**YEAR  
warranty



Dual /Y-block Flowmeter Assembly



Dual Outlet with Flowmeter

## Non-Removable Flowmeter


## Accessories:



Reusable Humidifier



Disposable Humidifier



Swivels/Christmas Trees



Power Takeoff

# Dial Style



Oxygen Dial Flowmeter



Medical Air Dial Flowmeter

## Features:

- Available in Oxygen in 25, 15, 8 and 4 LPM
- Available in Medical Air in 25 and 15 LPM
- Smallest and lightest in the industry
- Smallest increment is 1/32
- Color-coded labeling for specific gas
- Encased in an anodized aluminum housing
- Easy to read window

## Thorpe Style Technical Specifications

Physical			
Dimensions:	5.6" H x 2.7" W		
Weight:	6 oz		
Flow Range	Oxygen:	0 - 15 LPM	Medical Air:
		0 - 8 LPM	
		0 - 3 LPM	
		0 - 1 LPM	
		0 - 200 CC	
<ul style="list-style-type: none"> <li>• Back pressure compensated so flow accuracy is not affected by restrictions</li> </ul>		<ul style="list-style-type: none"> <li>• Calibrated at 50 psi</li> <li>• Serialized for ease of identification</li> </ul>	

## Dial Style Technical Specifications

Flow Range	Oxygen:	0 - 25 LPM	0.5, 1, 2, 3, 4, 5, 6, 8, 10, 15, 25
		0 - 15 LPM	0.5, 1, 1.5, 2, 3, 4, 5, 6, 8, 10, 15
0 - 8 LPM		0.5, 1, 1.5, 2, 2.5, 3, 4, 5, 6, 7, 8	
0 - 4 LPM		0, 0.03, 0.06, 0.12, 0.25, 0.5, 0.75, 1, 1.5, 2, 3, 4	
	Medical Air:	0 - 15 LPM	0.5, 1, 1.5, 2, 3, 4, 5, 6, 8, 10, 15
		0 - 25 LPM	0.5, 1, 2, 4, 6, 8, 10, 12, 15, 20, 25
<ul style="list-style-type: none"> <li>• Calibrated at 50 psi</li> </ul>		<ul style="list-style-type: none"> <li>• Serialized for ease of identification</li> </ul>	

## Model Number FMX-XXY-XX-Y

