



EMERGENCY OXYGEN UNITS & CPR ADMINISTRATION EQUIPMENT

[Emergency Oxygen Units](#)

[Cylinders](#)

[Regulators](#)

[Masks](#)

[Home](#)

[AHA Emergency Oxygen Guidelines](#)

[Contact LIFE®](#)

[About Us](#)

[Refill Information](#)

[Maintenance](#)

[Regulations](#)

[Warranty](#)

[Go To LIFE® StartSystem](#)

[Go To LIFE® SoftPac™](#)



LIFE is Time™
When minutes count
EMERGENCY OXYGEN
can be the difference between
LIFE OR DEATH.

- **Simple and EASY TO USE,**
External Two-step Instructions and Clear Cover
Inhalator: administer supplemental oxygen for a breathing victim
Resuscitator: oxygen enriched mouth-to-mask CPR for a non-breathing victim
- **Easily Portable** with Convenient Handle
20" x 6" x 9" - 16lbs - Aluminum Cylinder
- **Wall Mount** for Recognition and Accessibility, always available on display for easy access
- **Simple "On-Off" Control Lever**
- **Constant Reading Supply Gauge**
(other units must be turned on to read gauge)
- **90 Minute Supply @ 6LPM - 566 liters Oxygen**
- **Resuscitation LIFE® CPR Mask** with One-Way Valve fits both adult and child
- **Cylinder is easily disposable/replaceable** or refillable at local gas distributors.
- **Tight Seal Protects Unit from Dust and Moisture**
- **Units Come Ready-to-Use** full of oxygen with wall brackets and screws
- **5 Year Warranty**



This unit could save your life.



Make LIFE® OxygenPac part of your first aid and safety program.

[DEMO VIDEO](#)



LIFE® CPR Mask
invertably fits both ADULT and CHILD

IDEAL LOCATIONS FOR THE LIFE® OxygenPac
...WHEREVER THERE ARE PEOPLE!

For Employees & First Responders
INDUSTRIAL/FACILITY
PUBLIC FACILITIES
ENTERTAINMENT / RECREATION OFFICES



LIFE® OxygenPac is attractive. Prominently displayed, it makes the simple but strong statement,
"This Organization Cares For Its Employees And Visitors."

Uncompromisingly manufactured with quality materials, technology, and workmanship

MODEL	FEATURES	Regulators
LIFE-612	6&12 LPM - " NORM " & " HIGH " (2 settings) provides to recommended 100% inspired O2	
LIFE-025	0 to 25 LPM Variable Flow (only for trained EMTs) 12 flow settings - 0, .25, .5, 1, 2, 3, 4, 6, 8, 10, 15, 25	
LIFE-6FF	6 LPM (Liters Per Minute) Fixed-Flow	
LIFE-101	Replacement Cylinder Assembly complete with cylinder, valve, gauge & O2 fill	

[FDA Regulation](#) minimum capability is 6 LPM flow rate.
LIFE® 612 models deliver both the minimum of 6 LPM and the AHA recommended 100% inspired oxygen at 12 LPM, in just two simple settings which read "NORM" & "HIGH".
Only LIFE Corporation offers the simple 6&12 LPM flow rates.

American Heart Association Guidelines for CPR Oxygen

- ▶ "... give 100% inspired oxygen during Basic Life Support ... and Advanced Cardiac Life Support ... (BLS & ACLS) ... as soon as it becomes available."
- ▶ "... use supplemental oxygen...at a flow rate of 10-12 LPM ... to enable delivery of 100% inspired oxygen" (initially).