



## Keeler PSL Portable Slit Lamp-available with Camera options

10X & 16 X Mags, includes Charging Base

### PSL Portable Slit Lamps

Versatility and portability are key benefits of the Keeler PSL one and PSL Classic Slit Lamps. Keeler designed and manufactured with the strength to withstand the daily traumas of the demanding ophthalmic environment, the PSL is lightweight, durable and robust.

**Advanced optics:** 10x magnification, PSL one  
10x and 16x magnifications, PSL Classic

Top quality multi-layered glass is crafted into our Keeler multi-element and multi-layer coated optical system to give maximum performance. Clarity is guaranteed.

#### Controllable Illumination

Proven halogen illumination system. Illumination levels are controlled by a rheostat giving you control of your examination parameters. From zero to maximum, you choose.

#### Big Slit Lamp features, portable usability

Unique features such as the fixation targets, and the 1mm square light patch for assessing anterior chamber flare set the new standard for Portable Slit Lamps.



#### Carrying Case

Protect your investment with a rugged aluminum carrying case - perfect when your Slit Lamp is not in use or when you are travelling.



PSL one



PSL Classic

[www.keelerusa.com](http://www.keelerusa.com)

20



Portable-with Camera-iPhone adapter-Photos anywhere!

