

100 & 4200 SERIES



Pneumatic lift stool
Model 4246



Screw lift stool
Model 110

Pneumatic lift stool
Model 4240



Screw lift stool
Model 112

Rehance[®]

100 & 4200 SERIES ECONOMY STOOLS



Model 4240

Features a 15" diameter seat with a 2" high density polyfoam cushion. Seat height range is 18" to 26". Height is controlled with a convenient lever. Lift the lever and a pneumatic device assists in elevating the seat as you take your weight off the stool. To lower while seated, raise the lever and the stool descends gently to desired height.



Model 4246

A contoured, padded backrest adjusts horizontally and vertically for optimum back support. Seat is 15" in diameter with a 2" high density polyfoam cushion. Seat height range is 18" to 26". Height is controlled with a convenient lever. Lift the lever and a pneumatic device assists in elevating the seat as you take your weight off the stool. To lower while seated, raise the lever.



Model 110

Five-legged, 19" diameter base is standard with glides available as an option. Seat is 14" in diameter with a 2" polyfoam cushion. Height range is 18" to 26". Seat height is controlled by rotating the seat. To raise, hold the lift cap stationary. Rotating the seat clockwise raises the seat. To lower, hold the lift cap stationary and rotate the seat counterclockwise.



Model 112

A contoured padded backrest adjusts horizontally and vertically for optimum back support. Seat is 14" in diameter with a 2" polyfoam cushion. Height range is 18" to 26". Five-legged, 19" diameter base is standard with glides available as an option. Seat height is controlled by rotating the seat. To raise, hold the lift cap stationary. Rotating the seat clockwise raises the seat. To lower, hold the lift cap stationary and rotate the seat counterclockwise.

Swivel 360°

Seats on all 4200 and 100 Series Stools swivel 360° independently of the base, so height isn't affected.

Finished for Durability

All exposed metal, with the exception of the aluminum base casting, is chrome-plated steel. Base casting and legs are finished in tough, black epoxy enamel.

Upholstery

All stools come with black as standard, but you can match your office or treatment room decor with our selection of 14 easy-to-maintain vinyl colors. Other vinyl colors and nylon tweed inserts are available at additional cost. All materials are fire retardant. To request a color selection card, contact your dealer or call Reliance at 1-800-735-0357.

Stool Footring

A footring is available as an option. Order one with your stool or add at a later date. Approximate shipping weight: 2 lbs.

Dimensions and Specifications

Shipping weights: Model 4240 (23 lbs.), Model 4246 (26 lbs.), Model 110 (19 lbs.), Model 112 (26 lbs.).



Precision by Tradition